

EQUIPMENT USAGE

Cow Lifting Machine

Veterinarian often come across Cows and Buffalos with the “downer syndrome” unable to stand up needs to be treated by lifting the heavy animals to avoid complications of bed sores and nervous debility by improving circulation by massage. This Animal lifting Apparatus helps the Veterinaries to lift the recumbent animal for the treatment as well as close examination and evaluation.



Equipment to be taken to the site of recumbent animal and place it in such a way that the animal is in the middle of the equipment.

Fix the belt passing through the ventricle surface of the body from

Cell: 98780 - 72154

Email: sunnygenerators@gmail.com

Product Specification

Capacity	-	1 Ton
Usage / Application	-	Cow Lifting
Design Type	-	Trolley
Operating Height	-	2 Feet
Material	-	Mild Steel
Features	-	Rust proof
Country of origin	-	India

SUNNY ENGINEERING WORK





front portion to back ward one by one fixing the 8 hooks on either side and secure it well.

- ☞ Slowly lift the animal by operating the handle of the pulley rotating from left to right.
- ☞ Animal to be lifted to the height so that it can stand on its four limbs.
- ☞ Massage the body to increase the circulation.
- ☞ Veterinarian to decide course of treatment observing the improvement in bearing the load on itself.

MAINTAINENCE

For its long shelf life the equipment is to be properly maintained in the premises.

- ☞ The equipment is to be kept in a clean and dry condition.
- ☞ The belt is to be washed and dried after the usage and stored in a bag.
- ☞ Gear box system needs greasing once in a month.
- ☞ Trolley wheel to be greased regularly once in a month.